2026-06-11

Field Safety Notice

Date: 11 June 2026
To the intention of: Users and customers of Human Med AG
FSCA- Ref.: 2026-003

1. Details regarding the medical device concerned

Device name:	WAL irrigation and aspiration cannula
REF no.:	503001
Batch no.:	938823
Date of manufacture:	2025-10-27
Manufacturer:	Human Med AG

2. Reason for the Field Safety Notice

This issue was identified following a customer complaint. It was found that, for the afore-mentioned batch, the information on the packaging label does not match the marking on the cannula tube.

Specifically, the label on the packaging indicates the WAL irrigation and suction cannula REF 503001 with a diameter of 3.8 mm and 4 standard suction holes. The marking directly on the cannula tube, however, indicates the cannula variant REF 503002.

The WAL irrigation and suction cannula REF 503002 is a 'RAPID' variant; it has identical dimensions and the same hole geometry as the 'standard' variant REF 503001. Unlike the "standard" variant REF 503001, the suction holes of the "RAPID" variant REF 503002 are specifically designed for indications/treatments requiring more intensive suction.

Due to this discrepancy in the product labelling directly on the product and on the packaging label, clear identification of the cannula cannot be guaranteed, which may lead to misunderstandings regarding the functionality and use of the cannula.

3. Potential Risk

The use of a 'RAPID' cannula variant exposes the adipose tissue to higher shear forces during suction. If the user fails to recognize that the cannula he is using is a 'RAPID' variant and therefore does not adapt his technique accordingly, this may result in greater tissue trauma than intended. Consequently, increased tissue irritation or injury, haematomas and bleeding may occur.

4. Action required

- (1) Prior to use, check the information given on the packaging label of the above-mentioned cannula and on the cannula tube. If the REF numbers do not match, do not use the cannula!
- (2) Remove the incorrectly labeled cannula and mark it clearly to prevent further use.
- (3) If necessary, inform any internal departments or third parties to whom the affected cannula has been passed on.
- (4) Return the incorrectly labeled cannula to Human Med AG or to your supplier. For this, complete the product return form (Attachment 1) and enclose it with your return.

Please confirm receipt of this field safety notice and the return of the above-mentioned cannula using the enclosed confirmation form (Attachment 1) and send it to the email address below.

If you have any questions regarding this safety notice or the measures outlined, please contact our colleague:

Carina Buck

Tel.: 0049 (0)385 39570-25, Mobile: 0049 (0)151 57828292,

Email: buck@humanmed.com

Please contact your supplier for a refund of the original purchase price.

We sincerely apologize for any inconvenience caused and thank you for your understanding.

Kind regards

Bianca Tomuschat

Deputy PRRC (Art. 15 MDR)

**BITTE DIESES DOKUMENT AUSSEN AM VERSANDKARTON ANBRINGEN!
PLEASE FIX THIS DOCUMENT TO THE OUTSIDE OF THE TRANSPORT BOX!**

Bitte ausfüllen / Please fill out:

Kundeninformationen / Customer information

FSCA-Ref.: 2026-003

*Rücksendungsdatum / *Return date:

*Kundennummer / *Customer number:

*Adresse / *Address:

*PLZ, Stadt, Land / *Zip code, City, Country:

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***Date, Signature to confirm that the product has been cleaned and disinfected**

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Date, Signature of the recipient:

Dear Human Med Team,

I hereby confirm that I have read and understood your **Field Safety Notice dated 11 June 2026** relating to the

WAL irrigation and aspiration cannula, REF 503001, LOT 938823

Furthermore I confirm that the affected cannula(s) have been

Quantity: _____ pcs.

<input type="checkbox"/>	Disposed of	Date:
--------------------------	-------------	-------

<input type="checkbox"/>	Returned to the supplier	Date:
	Supplier's name and address	

<input type="checkbox"/>	Returned to Human Med AG	Date:
--------------------------	--------------------------	-------

Date:	
Name of the healthcare facility:	
Address:	
Name, first name:	
Signature:	