Urgent Field Safety Notice**Eckert & Ziegler
BEBIG GmbH**Robert-Rössle-Str. 10
13125 Berlin
Germany
www.bebig.com
SRN: DE-MF-000005760Contact: Bernd Schumacher
Tel: +49(0)30941084-130
E-Mail: customer.service@bebig.com

Subject: Wrong Certificate No.:3659625 for I25.S16

Commercial Name of Affected Product: I25.S16

UDI-DI 04049223100292

Reference: 2025OPH-2

Date of Notification: 09. Dec 2025

Type of Action: Replacement of certificate

Description of Problem:

Due to a human failure, the indicated total and apparent activity on certificate no.: 3659625 for lot number 0039845/10, order number BBD007981 were incorrect. This error concerns a single delivery. Using the incorrect value for apparent activity would result in an overdose of approx. 10%. An electronic version of the corrected certificate is attached to this information. A paper version of the corrected certificate will be provided via postage.

The recipient of this information is requested to replace the incorrect certificate and assure that only the corrected version is used for the planning of treatments.

This notice needs to be passed on to all those who need to be aware within your organisation or to any organisation where the devices have been transferred to.

Please reply to this email by 15. December 2025 confirming that you have received and understood this information and that you have forwarded it to the people which already have received the named products or will receive them.

The undersigned representative from Eckert & Ziegler BEBIG GmbH confirms that this notice has been notified to the appropriate Regulatory Agency.

We sincerely apologize for any inconvenience and thank you in advance for your cooperation. For further information please feel free to contact us.

Kind regards,

Bernd Schumacher
Head of Quality Management

Managing Directors:
Katrin Antonenko
Dr. Dirk W. BeckerBerlin Office:
Amtsgericht Charlottenburg
Reg.-Nr. HRB 42949 B
USt.-ID-Nr. DE 137169788Bank Account:
Commerzbank AG
IBAN DE29 1204 0000 0042 4648 00
BIC COBADEFFXXX

Urgent Field Safety Notice

Eckert & Ziegler
BEBIG GmbH

Robert-Rössle-Str. 10
13125 Berlin
Germany
www.bebig.com
SRN: DE-MF-000005760

CONFIRMATION

**Urgent Field Safety Notice from
Eckert & Ziegler BEBIG GmbH,
Reference 2025OPH-2**

This is to confirm that we have received and understood the Field Safety Notice. It was forwarded inside our clinic to the respective personnel.

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**EN Certificate for sealed radioactive sources
(IsoSeed®)**

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This certificate corrects the previous revision,
signed 26. nov. 2025

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Brachytherapy source

Radioactive sources **IsoSeed®** type **I25.S16** are intended for brachytherapy of solid tumours, especially in ophthalmology.

IsoSeeds are cylindrical, sealed sources. Their source type is specified by the product code. The seeds have dimensions of 4.5 mm in length and 0.8 mm in diameter. The encapsulation is a titanium tube with wall thickness of 0.05 mm, laser welded on both ends. Inside, the radionuclide I-125 is bound on a cylindrical carrier. The latter is working as X-ray marker additionally.

The half-life of I-125 is 59.41 days.

IsoSeeds are offered according to different source strength ranges. For that purpose, the reference air-kerma rate of every individual seed is measured. The reference date is referred to 12:00 noon MEZ/MESZ. The statement of "apparent activity" is based on dividing the measured quantity by constant A (see table).

product code	constant A in $\mu\text{Gy}\cdot\text{h}^{-1}\cdot\text{mCi}^{-1}$ (apparent activity)	constant B (contained activity)
I25.S16	1.27	1.3

Really contained activity results from "apparent activity" by multiplication with constant B (see table).

The measuring device used by Eckert & Ziegler BEBIG GmbH is periodically recalibrated. The calibration is traceable to the primary standard of the PTB, Germany.

The uncertainty of a single seed within a lot is at maximum $\pm 5\%$ ($k=2$).

Please use the following factor for the conversion of the indicated activity into MBq:

1 mCi = 37 MBq

Before shipment, all sources are checked for leak tightness, according to ISO 9978 by an immersion test in a suitable liquid at 70°C for 30 min.

They are provided non-sterile and have to be sterilized before application.

IsoSeeds are manufactured in accordance with a strict quality assurance program and in compliance with the standard DIN EN ISO 13485.

They are CE-marked according to 93/42/EEC.