

FSN Ref: FSCA_EGFR_MULTI_01_2025
 FSCA Ref: FSN_EGFR_MULTI_01_2025

Date: 21. 03. 2025

Urgent Field Safety Notice **gb ONCO EGFR**

For Attention of*: Customers and distributors.

Contact details of local representative	
Name of local representative:	Jitka Kašparová
e-mail:	jitka.kasparova@generi-biotech.com
phone:	+420 734 334 449
Manufacturer address:	GENERI BIOTECH s.r.o. U Fotochemy 1763 500 02 Hradec Králové – Pražské Předměstí

A. Information on Affected Devices*	
1.	Device Type(s)*
	IVD kit, non-sterile
2.	Commercial name(s)
	gb ONCO EGFR
3.	Unique Device Identifier(s) (UDI-DI)
	N/A
4.	Primary clinical purpose of device(s)*
	gb ONCO EGFR diagnostic kit enables detection of EGFR mutations in exons 18, 19, 20 and 21 in cancer patients.
5.	Device Model/Catalogue/part number(s)*
	3282-024
6.	Software version
	N/A
7.	Affected serial or lot number range
8.	200104005
9.	Associated devices
	N/A

B. Reason for Field Safety Corrective Action (FSQA)*	
1.	Description of the product problem* When individual reactions are prepared, patient outcomes cannot be evaluated, the analysis is non-valid.
2.	Hazard giving rise to the FSQA* Invalid analysis, delay in results release
3.	Probability of problem arising The probability is very high. Invalid analysis only occurs if the kit is used to prepare individual reactions (for each patient). Once the kit is used once for all 24 reactions, the analysis is valid and correct results are obtained.
4.	Predicted risk to patient/users Delay in results release.
5.	Further information to help characterise the problem N/A
6.	Background on Issue Customer complaint.
7.	Other information relevant to FSQA N/A

C.					
D. Type of Action to mitigate the risk*					
1.	Action To Be Taken by the User* <p> <input checked="" type="checkbox"/> Identify Device <input type="checkbox"/> Quarantine Device <input checked="" type="checkbox"/> Return Device <input checked="" type="checkbox"/> Destroy Device </p> <p> <input type="checkbox"/> On-site device modification/inspection </p> <p> <input type="checkbox"/> Follow patient management recommendations </p> <p> <input checked="" type="checkbox"/> Take note of amendment/reinforcement of Instructions For Use (IFU) </p> <p> <input type="checkbox"/> Other <input type="checkbox"/> None </p> <p> <i>Provide further details of the action(s) identified.</i> Do not use the kit to prepare single reactions. </p>				
2.	<table border="1"> <tr> <td>By when should the action be completed?</td> <td>Immediately.</td> </tr> </table>	By when should the action be completed?	Immediately.		
By when should the action be completed?	Immediately.				
3.	<table border="1"> <tr> <td>Particular considerations for:</td> <td>IVD</td> </tr> <tr> <td> Is follow-up of patients or review of patients' previous results recommended? <i>Provide further details of patient-level follow-up if required or a justification why none is required</i> </td> <td>No</td> </tr> </table>	Particular considerations for:	IVD	Is follow-up of patients or review of patients' previous results recommended? <i>Provide further details of patient-level follow-up if required or a justification why none is required</i>	No
Particular considerations for:	IVD				
Is follow-up of patients or review of patients' previous results recommended? <i>Provide further details of patient-level follow-up if required or a justification why none is required</i>	No				

FSN Ref: FSCA_EGFR_MULTI_01_2025
FSCA Ref: FSN_EGFR_MULTI_01_2025

4.	Is customer Reply Required? * (If yes, form attached specifying deadline for return)	No
5.	Action Being Taken by the Manufacturer <input checked="" type="checkbox"/> Product Removal <input type="checkbox"/> On-site device modification/inspection <input type="checkbox"/> Software upgrade <input checked="" type="checkbox"/> IFU or labelling change <input type="checkbox"/> Other <input type="checkbox"/> None <i>Provide further details of the action(s) identified.</i>	
6.	By when should the action be completed?	31. 03. 2025
7.	Is the FSN required to be communicated to the patient /lay user?	No
	If yes, has manufacturer provided additional information suitable for the patient/lay user in a patient/lay or non-professional user information letter/sheet?	Not assigned for this FSN

E. General Information*		
1.	FSN Type*	New
2.	For updated FSN, reference number and date of previous FSN	N/A
3.	For Updated FSN, key new information as follows	N/A
4.	Further advice or information already expected in follow-up FSN? *	No
5.	If follow-up FSN expected, what is the further advice expected to relate to:	N/A
6.	Anticipated timescale for follow-up FSN	N/A
7.	Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	GENERI BIOTECH s.r.o.
	b. Address	U Fotochemy 1763 500 02 Hradec Králové – Pražské Předměstí
	c. Website address	www.generi-biotech.com
8.	The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *	
9.	List of attachments/appendices:	N/A
10.	Name/Signature	Mgr. Jitka Kašparová, Ph.D. <small>Digitálně podepsal Jitka Kašparová Datum: 2025.03.21 00:06:34 +01'00'</small>

FSN Ref: FSCA_EGFR_MULTI_01_2025
FSCA Ref: FSN_EGFR_MULTI_01_2025

Transmission of this Field Safety Notice

<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p>
--

<p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p>
--

<p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p>
--

<p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.*</p>
--

Note: Fields indicated by * are considered necessary for all FSNs. Others are optional.