

MIKROGEN GmbH | Anna-Sigmund-Str. 10 | 82061 Neuried | Germany

Contact
Martina Wild
vigilance@mikrogen.de
Phone +49 89 54801-145
Fax +49 89 54801-100

date
Neuried, 19.03.2026

page
1 von 3

Urgent safety information
(reply required, page 3)
due to possible false-negative results for influenza A virus

Products:

***ampliCube* Respiratory Viral Panel 5 (article-no. 50152)**
***ampliCube* Respiratory Flu & SARS-CoV-2 (article-no. 50163, 50164)**
***ampliCube* Respiratory Flu & SARS-CoV-2 LC (article-no. 50173, 50174)**

Subject: Mutations in the Influenza A virus (H1N1) pdm09 subtype may lead to false-negative results.

Dear Partner,

during the current 2025/26 season of acute respiratory infections, an influenza A virus (H1N1) pdm09 variant is circulating in the field carrying mutations in the target region used for the detection of influenza A virus. These mutations may affect the binding affinity of the primers used in the influenza A virus channel. Depending on the sampling and extraction efficiency, this may lead to a false-negative result in patient samples with a low viral load.

We would like to emphasise that our product *ampliCube* Respiratory Viral Panel 1 (article-no. 50102) is not affected by this issue, as the mutated variant of the subtype continues to be detected via the H1N1 channel.

Recommendations for you as a partner and your customers:

- Based on current knowledge, you and your customers can assume that no action is required with regards to historical data. Test results must always be considered in conjunction with the clinical picture. In recent cases, retesting may be considered for high-risk patients, e.g. those at risk of Reye's syndrome, at the discretion of the treating physician.

- Please inform us of your stock as well as your customers stock of the above-mentioned products. We will replace these accordingly:
 - o *ampliCube* Respiratory Viral Panel 5 with *ampliCube* Respiratory Viral Panel 1 (article-no. 50102) and *ampliCube* Respiratory Viral Panel 3 (article-no. 50122).
 - o *ampliCube* Respiratory Flu & SARS-CoV-2 (LC) for *ampliCube* Respiratory Viral Panel 1 (article-no. 50102) and *ampliCube* Coronavirus SARS-CoV-2 (article-no. 50143)
- Please destroy any kits of the above-mentioned products currently in your stock. Please ensure that your customers also destroy all affected kits in stock.

Please ensure that all customers working with the above-mentioned product are made aware of this safety information. As part of this measure, the products will be updated and made available as soon as possible.

The BfArM has been informed of this measure.

We ask you to confirm receipt of this safety information by completing the reply form on page 3. You can send this reply to us by email to vigilance@mikrogen.de or by fax to +49-89-54801-100. We kindly ask for your answer by 30th April 2026.

We sincerely apologise for any inconvenience caused and thank you for your support.

Should you have any questions regarding this measure, please do not hesitate to contact us.

Kind regards,

Vigilance

Response from the distribution partner regarding the safety information for

ampliCube Respiratory Viral Panel 5 (article-no. 50152)
ampliCube Respiratory Flu & SARS-CoV-2 (article-no. 50163, 50164)
ampliCube Respiratory Flu & SARS-CoV-2 LC (article-no. 50173, 50174)

We kindly ask you to send us your response by email to vigilance@mikrogen.de or by fax to +49-89-54801-100. We would appreciate your response by 30 April 2026. Thank you very much!

| | |
|--|--|
| 1. I have read and taken note of the recommended measures. | <input type="checkbox"/> yes <input type="checkbox"/> no |
| 2. I have forwarded the recommended measures to my affected customers and received feedback regarding kit numbers. | <input type="checkbox"/> yes <input type="checkbox"/> no |
| 2. Number of kits (in your stock and in your customers' stock): | |
| <i>ampliCube Respiratory Viral Panel 5</i> | |
| <i>ampliCube Respiratory Flu & SARS-CoV-2</i> | |
| <i>ampliCube Respiratory Flu & SARS-CoV-2 LC</i> | |
| 3. I confirm that the kits listed above have been properly disposed of. | <input type="checkbox"/> yes <input type="checkbox"/> no |

We will contact you with further details regarding the dispatch of the replacement kits.

If no, or if you do not wish an exchange, please explain why:

Please specify the countries to which the affected kits were distributed:

Name _____

Position _____

Company _____

Street/Number _____

Postcode/location _____

Date and signature _____