

For the attention of
Distributors and Users

11. Juni 2026 / BA
esther_bauer@fehling-instruments.de
www.fehling-instruments.de

FS-26-06

Urgent safety information

RECALL

regarding

MOA-2 Biopsy forceps, single use, sterile, 1,6 x 800 mm

Identification of affected medical devices:

REF: MOA-2
LOT: OG2
Description: Biopsy forceps, single use, sterile, 1,6 x 800 mm
Basic-UDI: 4058749-FI10-BQ
UDI-DI: 4058749000023

Description of the problem including the identified cause:

The biopsy forceps MOA-2, LOT # OG2, were manufactured out of specifications. The length of the biopsy forceps is only 510 mm instead of 800 mm as indicated on the label.

If these products are used, the biopsy may not be able to be performed as planned due to the length of 510 mm of the instrument.

What measures are to be taken by the addressee?

Please return all concerned products to the manufacturer via your local distributor.

Measures taken by FEHLING INSTRUMENTS GmbH

FEHLING INSTRUMENTS will replace the biopsy forceps by biopsy forceps with correct length.

Disclosure of the information described herein:

Please ensure in your organization that all users of the above-mentioned products and other persons to be informed are made aware of this Urgent Safety Information.

If you have given the products to a third party, please forward a copy of this information or inform the contact person indicated below.

Please keep this information at least until the action has been completed.

The Federal Institute for Drugs and Medical Devices has received a copy of this "Urgent Safety Information".

Contact:

If you have any questions about the process, please contact the person below:

FEHLING INSTRUMENTS GmbH
Seligenstädter Str. 100
63791 Karlstein, Germany
www.fehling-instruments.de

Esther Bauer
Tel.: +49 6188-95740
vigilance@fehling-instruments.de

Kind regards

FEHLING INSTRUMENTS GmbH

Esther Bauer
Head of Quality Management/ Regulatory Affairs

CONFIRMATION FORM

Please return by e-mail to

vigilance@fehling-instruments.de

FSN NO.:	FS-26-06
FSN Date:	2026-06-10
Product	MOA-2 Biopsy forceps, single use, sterile, 1,6 x 800 mm

CUSTOMER

Name	
Address	
Contact	
Phone number	
Email	

Measures to be carried out by the customer

- I acknowledge receipt of the FSN and confirm that I have read and understood its contents.
- The information and required actions were brought to the attention of all relevant users and implemented.
- The concerned products will be returned to the distributor.

Name	
Date	
Signature	