

Ralf Quentin
Department Manager Quality Assurance
& Medical Device Surveillance

Tel.: +49 241 / 91 30 0
E-Mail: rquentin@vygon.com

Aachen, 17.03.2026

Urgent FSMA on lifecath CT PICC easy – Recall 011

As legal manufacturer of the affected article, we would like to inform you that we have identified a batch-related manufacturing problem regarding the size of the catheter tube.

Due to a manufacturing fault, in some cases a 5F-catheter was placed in the 4F-catheter set. This catheter is too big and cannot be inserted through the included 4F-sheath.

Product Description	Code No.	Batch No.
lifecath CT PICC easy	V021292414	171125GE

We therefore ask you to immediately withdraw this batch from services and to return affected products according to our instructions below.

Please acknowledge receipt of this letter and complete and return the attached form, indicating the quantities withdrawn from your institution.

Your National Competent Authority has been informed about this Recall.

If you should require further information, please contact your local VYGON distributor.

We apologize for any inconvenience this recall may cause.

Yours sincerely,

Signed by Ralf Quentin

Department Manager Quality Assurance & Medical Device Surveillance

Urgent FSCA on lifecath CT PICC easy – Recall 011

ACKNOWLEDGMENT AND CUSTOMER RESPONSE FORM

Please complete and return this form by e-mail to vigilance@vygon.com

Name and address of the institution:	
Full name of the person to contact:	
Function:	
☎ Phone number:	
✉ E-mail:	

We acknowledge receipt of the above FSCA and that the information contained in this field safety notice has been shared with all recipients / end users of above-mentioned products within our organization.

Please tick the appropriate box:

We have the following products available on stock:

Code No.	Batch No.	Number of units used	Number of units blocked
V021292414	171125GE		

Signature and Date: _____