



LTA MEDICAL  
27 Rue Albert Einstein,  
93000 Bobigny  
01.48.58.33.55

## Field Safety Notice

Single use endoscope cleaning brush  
Recall

January 2025

For Attention of: Healthcare professional, Vigilance correspondent

The purpose of this letter is to inform you that LTA MEDICAL is conducting a recall of its brushes.

### **Intended use:**

The brushes are single use and are intended for cleaning flexible endoscopes.

### **Problem Description:**

This recall follows reports of the brush detaching from its stem.

### **Risks for the patients:**

The use of these brushes may lead to:

- The brush getting stuck inside the endoscope.
- The risk of the brush becoming lodged in the patient's body during the use of the endoscope, potentially leading to infections.

### **Affected Products:**

Below is the list of relevant references.

Reference	Lot number	GTIN
100-1824/50	24300356	3760274483875
100-1824/50S	24490097	3760274488245
200-1824/5	24370036	3760274487286
200-1824/50	24250255	3760274483967
	24260007	
	24270186	
	24370015	
200-1824/50	24430055	3760274483967
BSC-1824	24470184	3760274485220
	24480378	
	24490044	
BRONC-STERIL	24500091	3760274487583
D1111825ND/25	24390035	3701571301493
D111825ND/25	24380158	3760274484865
	24420142	
D11-2632/25	24430143	3760274484377
	24440144	
D11-3748/25	24420037	3760274484391
D-1824/25	24270165	3760274484346
	24370019	
	24380294	
D-2632/25	24420160	3760274484360
D-3748/25	24450027	3760274484384
DI-3748/25	24420049	3760274484483
EB100-18/25	24420038	3760274484780
	24430045	
	24470223	
EB-1824/25	24290027	3760274484117
	24340172	
	24420183	
EBI-1824/25	24370021	3760274484148
	24450193	
	24470183	
	24480008	
NGEE-1824	24510048	3760274485299
	24400066	
	24410166	
R-1824/25	24340170	3760274484179
	24350099	
RI-1824/25	24290032	3760274484209
	24400219	
	24490283	

SD11-3748/25	24440165	3760274484575
SD-1824/25	24250254	3760274484520
	24420170	
	24360100	
SD-3748/25	24450030	3760274484568
SEB100-18/25	24340137	3760274484827
SEB-1824/25	24290034	3760274484230
	24350102	
	24420185	
SEBI10018/25	24310104	3760274484834
	24360018	
SEBI10018/25	24430134	3760274484834
	24480027	
SR-1824/25	24310118	3760274484292
	24380300	
	24420186	

Please verify whether you have any of the products mentioned.

**Instructions for the Customer:**

**Identification and Quarantine:** As soon as possible, identify the products in your inventory and quarantine the affected ones.

**Return of Products:** Return the affected products to LTA MEDICAL by sharing the information to our QARA team through mail address [regulatory.lta@getinge.com](mailto:regulatory.lta@getinge.com) . Our back office will organize the replacement of the products.

**Acknowledgment:** Please complete and return the attached acknowledgment form, even if you do not have any affected products in your inventory.

**LTA MEDICAL has notified the Competent Authorities of this action.**

We sincerely apologize for any inconvenience this may cause. Patient safety remains our top priority, and we thank you in advance for your prompt action. For any questions related to this communication, please contact: [regulatory.lta@getinge.com](mailto:regulatory.lta@getinge.com)

Sincerely,

Jerome POLENNE

Managing Director / PRRC

**Attachments:**

Acknowledgment form



LTA MEDICAL  
27 Rue Albert Einstein,  
93000 Bobigny  
01.48.58.33.55

### Acknowledgment of receipt

I, the undersigned, hereby acknowledge receipt of the urgent safety notification and confirm that I have fully understood the actions required. Furthermore, I certify that this safety notification has been duly distributed to the relevant personnel.

Reference	Lot number	Quantity

Country : .....

Facility : .....

Date : .....

Name (function) : .....

Signature : .....

Document to be returned to:

[regulatory.lta@getinge.com](mailto:regulatory.lta@getinge.com)