

Aachen, 04.07.2025

### **Urgent FSCA on lifecath CT PICC easy – Recall 010**

As legal manufacturer of the affected article, we would like to inform you that we have identified a batch-related manufacturing problem regarding the screw nut for connection between extension line and catheter tube.

Due to a manufacturing fault during the assembly of the extension line there is the possibility that a wrong size of screw nut was used, which may result in insufficient holding strength when connecting the catheter tube.

Product Description	Code No.	Batch No.
lifecath CT PICC easy	V021292214	230125GF, 210225GG
	V021292414	310125GD
	V021292514	260225GD

We therefore ask you to immediately withdraw these batches from services and to return them according to our instructions below.

Please acknowledge receipt of this letter and complete and return the attached form, indicating the quantities withdrawn from your institution.

Your National Competent Authority has been informed about this Recall.

If you should require further information, please contact your local VYGON distributor.

We apologize for any inconvenience this recall may cause.

Yours sincerely,



Signed by [REDACTED]

Department Manager Quality Assurance & Medical Device Surveillance

VYGON Germany GmbH  
Tel.: +49 (241) 91 30 0  
Fax: +49 (241) 91 30 106  
E-Mail: info-de@vygon.com  
WEEE-Reg-Nr.: DE49290163

Dr. Raymund Heiliger  
Geschäftsführer / Managing Director  
Ust-Ident-Nummer: DE 121 68 64 46  
Amtsgericht Aachen HRB 28087

Sparkasse Aachen  
Kto.-Nr. 53 884 - (BLZ 390 500 00)  
IBAN-Nr.: DE89 3905 0000 0000 0538 84  
BIC Code: AACSDE 33

Postbank Köln  
Kto.-Nr. 0211 680 508 - (BLZ 370 100 50)  
IBAN-Nr.: DE59 3701 0050 0211 6805 08  
BIC Code: PBNKDEFF370

## Urgent FSCA on lifecath CT PICC easy – Recall 010

### ACKNOWLEDGMENT AND CUSTOMER RESPONSE FORM

Please complete and return this form by e-mail to [vigilance@vygon.com](mailto:vigilance@vygon.com)

<b>Name and address of the institution:</b>	
<b>Full name of the person to contact:</b>	
<b>Function:</b>	
<b>☎ Phone number:</b>	
<b>✉ E-mail:</b>	

**We acknowledge receipt of the above FSCA and that the information contained in this field safety notice has been shared with all recipients / end users of above-mentioned products within our organization.**

**Please tick the appropriate box:**

**We have the following products available on stock:**

Code No.	Batch No.	Number of units removed
V021292214	230125GF	
	210225GG	
V021292414	310125GD	
V021292514	260225GD	

Signature and Date: \_\_\_\_\_