

Drägerwerk AG & Co. KGaA, 23542 Lübeck

To our customers of the  
Dräger Anesthesia Workstations Atlan

July 2026

## Second Addendum Important Safety Notice!

Possible failure of automatic ventilation  
due to a fault in the motor electronics

**Incorrect manufacturing date string within UDIs provided with our  
Atlan - Addendum Important Safety Notice!, dated April 2026**

Dear Sir or Madam

Due to a technical issue, the UDIs we provided with our last Field Safety Notice contain incorrect string section (11) for the manufacturing date.

| UDI                                                               |
|-------------------------------------------------------------------|
| (01)04048675556343(11) <b>240109</b> (21) ASSN-0311(93)8621300-15 |

Pic. 1 – incorrect string section (11)

To identify the affected devices properly, you may simply compare the serial number string section (21) with your serial numbers on the labelling of your devices (refer to Pic. 3).

| UDI                                                              |
|------------------------------------------------------------------|
| (01)04048675556343(11)240109(21) <b>ASSN-0311</b> (93)8621300-15 |


Pic. 2 – serial number string section (21)

Drägerwerk AG & Co. KGaA  
Moislinger Allee 53-55  
23558 Lübeck, Germany  
Postal address:  
23542 Lübeck, Germany  
Tel +49 451 882-0  
Fax +49 451 882-2080  
info@draeger.com  
www.draeger.com  
VAT no. DE135082211

Bank details:  
Commerzbank AG, Lübeck  
IBAN: DE95 2304 0022 0014 6795 00  
Swift-Code: COBA DE FF 230  
Sparkasse zu Lübeck  
IBAN: DE15 2305 0101 0001 0711 17  
Swift-Code: NOLADE21SPL

Registered office: Lübeck  
Commercial register:  
Local court Lübeck HRB 7903 HL  
General partner: Drägerwerk Verwaltungs AG  
Registered office: Lübeck  
Commercial register:  
Local court Lübeck HRB 7395 HL

Chairman of the Supervisory Board  
for Drägerwerk AG & Co. KGaA  
and Drägerwerk Verwaltungs AG:  
Stefan Lauer  
Executive Board:  
Stefan Dräger (chairman)  
Rainer Klug  
Gert-Hartwig Lescow  
Dr. Reiner Piske  
Anton Schrofner



Pic. 3 – Device label with serial number [SN]

We apologize for the inconvenience resulting from this error.

Makoto Onda

Global Product Manager, Atlan family  
Business Unit Therapy, Anesthesiology